

ALIGNER PLIERS



markhor™
cutters and pliers

Hole-Punch

6mm, round



Plier for punching out half-moon recesses in the area of brackets or bonded buttons. Also for removing protruding plastic edges to avoid damage to the soft tissue.



ORDER NO. OC 1001 A

Punch drop shape



Plier to add a drop shaped notch at the gingival margin. The drop shape enables the secure attachment of elastics - even for the patient.



ORDER NO. OC 1002 A

Horizontal Plier



Special plier for creating horizontal recesses for individual root torque. Also allows to improve the retention of aligners or retainers.



ORDER NO. OC 1003 A

Vertical Plier



Special pliers for over-correction of a rotation. Furthermore the vertical plier can eliminate in-and-out discrepancies.



ORDER NO. OC 1004 A

Point elevation plier

1,5mm



Special plier for creating punctiform elevations in the aligner. Allows a selective activation, which increases the pressure on the individual teeth.



ORDER NO. OC 1005 A

Reverse Plier



Special plier for inverting elevations in the aligner.



ORDER NO. OC 1006 A